

# USING NON-INVASIVE BRAIN STIMULATION TO MODULATE SUPRASPINAL CIRCUITS AFTER SCI

Jennifer Iddings, PhD & Anastasia Zarkou, PT, PhD

Postdoctoral Research Fellows

Crawford Research Institute

SCI Research Laboratory



Shepherd Center

2020 Peachtree Road, NW, Atlanta, GA 30309-1465

404-352-2020 [shepherd.org](http://shepherd.org)

# Non-Invasive Brain Stimulation (NIBS)

Two primary classifications of NIBS modalities:

1. Transcranial magnetic stimulation (TMS)

2. Transcranial current stimulation (tCS)



# Non-Invasive Brain Stimulation (NIBS)

Two primary classifications of NIBS modalities:

## 1. Transcranial magnetic stimulation (TMS)


- Neurophysiology measure: single & paired pulse TMS
- Neuromodulation tool: repetitive TMS (rTMS), paired associative stimulation (PAS)

## 2. Transcranial current stimulation (tCS)



# Basic Principles of TMS

- Electrical current in coil generates magnetic field
- Magnetic field passes through the skull unimpeded
- Alternating magnetic fields induce electrical current in brain




*George, Sci Am, 2003*



# TMS as Neurophysiology Measure

- Used to assess cortical & spinal circuitry
- Suprathreshold stimulation of the motor cortex elicits muscle response - motor evoked potential (MEP)
- Parameters of interest:
  - MEP amplitude
  - MEP latency
  - Silent period
  - Motor threshold (resting and active)



# Non-Invasive Brain Stimulation (NIBS)

Two primary classifications of NIBS modalities:

## 1. Transcranial magnetic stimulation (TMS)

- Neurophysiology measure: single & paired pulse TMS
- Neuromodulation tool: repetitive TMS (rTMS), paired associative stimulation (PAS)


## 2. Transcranial current stimulation (tCS)

- direct (tDCS)
- pulsed (tPCS)
- alternating (tACS)
- random noise (tRNS)



# Basic Principles of tCS

- **Subthreshold** stimulation modulating excitability of large brain regions
  - Does not elicit neuronal depolarization
- Mechanism: modulation of neuronal resting excitability
- Two categories:
  - Polarity dependent
    - Excitatory: anodal
    - Inhibitory: cathodal
  - Polarity independent
    - Excitation/inhibition are frequency and intensity dependent



*Adapted from Villamar et al., J Vis Exp, 2013*



# Basic Principles of tCS

## After-Effect Modulation

**Polarity  
Dependent**

**Polarity  
Independent**

---

*Waveforms adapted from: Jaberzadeh & Zoghi, Basic Clin Neurosci, 2013; Jaberzadeh et al., PLOS One, 2015*





# Application Basics of tCS

## tDCS Devices



## tCS Device tDCS, tPCS, tACS, tRNS



# Application Basics of tCS

- **Stimulation Parameters:**

- ✓ Current type: DC or AC
- ✓ Intensity: < 4mA
- ✓ Frequency: 0- 640Hz
- ✓ Duration: < 60 min/ day
- ✓ Ramp up/ down

Antal et al., Clin Neurophysiol, 2017; DaSilva et al., J Vis Exp, 2011; Woods et al., Clin Neurophysiol, 2016




# Application Basics of tCS

- **Stimulation Parameters**

- **Adverse Effects:**

- ✓ Itching
- ✓ Tingling
- ✓ Headache
- ✓ Discomfort
- ✓ Burning Sensations
- ✓ Skin redness
- ✓ Mild fatigue
- ✓ Skin burns



(Adapted from Palm et al., Brain Stimul, 2008)

., Brain Stimul, 2011)

Antal et al., Clin Neurophysiol, 2017; DaSilva et al., J Vis Exp, 2011; Woods et al., Clin Neurophysiol, 2016



# Application Basics of tCS


- Stimulation Parameters
- Adverse Effects
- **Contraindications:**
  - ✓ Implanted metallic device in the head
  - ✓ Uncontrolled seizures
  - ✓ Severe headache history
  - ✓ Scalp Lesions
  - ✓ Pregnancy
  - ✓ Pacemaker

Antal et al., Clin Neurophysiol, 2017; DaSilva et al., J Vis Exp, 2011; Woods et al., Clin Neurophysiol, 2016



# Application Basics of tCS:

## *Using Montages to Direct Effects*




<http://what-when-how.com/neuroscience/the-upper-motor-neurons-motor-systems-part-1/>



# Application Basics of TCS:

## *Montages for Upper Extremity Neuromodulation*




Adapted from [www.clinicalgate.com/epilepsy-8](http://www.clinicalgate.com/epilepsy-8)



# Application Basics of TCS:

## *Montages for Upper Extremity Neuromodulation*



**Both UEs**


*Adapted from [www.clinicalgate.com/epilepsy-8](http://www.clinicalgate.com/epilepsy-8)*





# Application Basics of TCS:

## *Montage for Lower Extremity Neuromodulation*



Adapted from [www.clinicalgate.com/epilepsy-8](http://www.clinicalgate.com/epilepsy-8)





# Advantages of tCS

tCS	
Focality	• Broad stimulation area
Cost	✓ More affordable
Portability	✓ Small footprint, easy to transport
Ease of Use	✓ Easily applied
Application	✓ Participant can move during stimulation




# Non-Invasive Brain Stimulation in Persons with SCI




# Maladaptive Cortical Plasticity Post-SCI:

## Reduced Motor Evoked Potential (MEP) Amplitude

1. Disruption of communication within the spinal cord
2. Decrease in cortical excitability





# Maladaptive Cortical Plasticity Post-SCI: Reduced Motor Evoked Potential (MEP) Amplitude



*Davey et al., Exp Brain Res, 1999*



# Reversing Maladaptive Cortical Plasticity of Thenar Muscles in Persons with SCI



*Hoffman & Field-Fote, Top Spinal Cord Rehabil, 2013*



# NIBS Studies in SCI Research: Recruitment Completed



# Conditioning Neural Circuits to Improve Upper Extremity Function


## Purpose:

- To examine how well people with tetraplegia are able to use their arms and hands after receiving fine motor training (FMT) combined with stimulation to increase brain excitability either directly or indirectly



Direct brain stimulation:  
Transcranial direct  
current stimulation


Indirect brain stimulation:  
Peripheral nerve  
somatosensory stimulation



# Conditioning Neural Circuits to Improve Upper Extremity Function

## Design

- Randomized controlled trial: each participant is randomly assigned to one of three groups: tDCS + FMT, PNSS + FMT, or sham tDCS + FMT





# Conditioning Neural Circuits to Improve Upper Extremity Function

## Outcome Measures:



- UE motor impairment and function
- UE strength and sensation
- Participant perception of function/participation
- Participant perception of quality of life
- Corticospinal excitability
- Spinal reflex excitability




# Enhancing Corticospinal Excitability to Improve Functional Recovery

## Purpose:


- To compare four different stimulation conditions for increasing brain excitability to determine which is best for helping people with tetraplegia improve their ability to use their arms and hands



### unihemispheric



### bihemispheric



# Enhancing Corticospinal Excitability to Improve Functional Recovery

## Design:

- Randomized crossover trial
- Each subject participates in a single session of each stimulation type

## Outcome Measures:

- Motor Control – timed tapping task
- Dexterity – nine hole peg test (9HPT)
- Strength – key pinch force
- Corticospinal excitability – motor evoked potentials




# NIBS Studies in SCI Research: Current & Upcoming



# Random Noise to Enhance Corticomotor Drive for Improved Hand Function


## Purpose:

- To compare the effects on cortical excitability, motor and sensory function of 3-day combined tRNS + fine motor training protocol to tDCS + fine motor training and to sham-stimulation + fine motor training protocols in individuals with tetraplegia.



# Random Noise to Enhance Corticomotor Drive for Improved Hand Function

Design:



Outcome Measures:

Cortical Excitability- Motor Evoked Potentials

Unimanual Function- Grasp and Release Test

Bimanual Function- CAHAI\_9

Sensory Function- Sensation subtest of GRASSP &  
revised Nottingham Sensory Assessment


Strength- Pinch and Grasp strength



# Enhancing Corticospinal Activation for Improved Walking Function

## Purpose:

- Determine if moderate-intensity, circuit training can improve walking-related outcomes
- Determine if the addition of non-invasive brain stimulation will result in greater improvements in function compared to circuit training alone



# Enhancing Corticospinal Activation for Improved Walking Function


## Design:

- Randomized controlled trial
- Each participant is randomly assigned to one of two groups: tDCS + training OR sham tDCS + training

## Outcome Measures:

- Walking Function
- LE strength
- Balance
- Spasticity
- Functional Strength

### CIRCUIT TRAINING



**Note:** Task modification and balance support provided as needed





# Acknowledgements

## Edelle Field-Fote, PT, PhD

Allison McIntyre, MOT, OT/L, CCRP

Barry McKay, BS

Beth Sasso, PT, DPT, NCS

Brandon Poe, MPA, CCRP

Cathy Furbish, PT

Cazmon Suri, MS

Evan Sandler, PT, DPT

Jasmine Hope, BS

Kyle Condon, PT, DPT

Marissa Mirecki, MOT, OTR/L

Nick Evans, MHS, ACSM, CEP

Ryan Koter, PT, DPT

Sarah Callahan, MOT, OTR/L

Stephen Estes, PhD

## Hulse Spinal Cord Injury Lab



*Eunice Kennedy Shriver* National Institute  
of Child Health and Human Development



Shepherd Center

[shepherd.org](http://shepherd.org)

# Questions?

